

iVascular[®]
therapies for living

navitian

Coronary microcatheter



The navigation you need



www.ivascular.global

Navitian, the needed microcatheter to arrive to the most challenging lesions.

Its exceptional design brings benefits as:

Outstanding trackability
through the small and
tortuous arteries due to
hydrophilic coating
HYDRAX PLUS


**Optimal flexibility,
pushability and resistance
to kinking** provided by a
full length variable braided
pattern

**Complete
guidewire control**
ensured by its
conical transition

**Superior crossability and
penetration capacity** due to
its optimized tip profile and
shape with double c-cut
radiopaque markers

navitian features

- > Over the wire catheter (OTW)
- > Guidewire compatibility 0.014"
- > Usable catheter length 135 and 150 cm
- > Radiopaque markers
- > Internal PTFE layer
- > Hydrophilic coating: 70 and 90 cm



Product with CE mark, certified by Notified Body 0318

References

Usable length 135 cm	MCC C14 135 001
Usable length 150 cm	MCC C14 150 001

The availability of each reference for the sale is linked to the authorization of commercialization in the country of destination

Distributed by:



Manufactured by:

Life Vascular Devices Biotech S.L
www.ivascular.global
info@ivascular.global