

















Presentation of the most popular ECG disposable electrodes of FARUM sp. z o.o.

<i>REF. NO.</i>	<i>IMAGE</i>	<i>DESCRIPTION</i>
FES-5541C		Adult, Holter / Stress, metal snap, Ag/AgCl sensor, solid gel, PE foam, size 55mm/41 mm. Packing: pouch – 50pcs., carton 2000 pcs.
FES-5136C		Adult, Stress/monitoring, metal snap, Ag /AgCl sensor, solid gel, PE foam, size 51mm/36mm. Packing: pouch – 50pcs., carton 2000 pcs.
FES-5036D		Adult, monitoring, metal snap, Ag /AgCl sensor, solid gel, PE foam, size 50mm/36mm. Packing: pouch – 50pcs., carton 2000 pcs.
FES-4543B		Adult, monitoring, metal snap, Ag /AgCl sensor, solid gel, PE foam, size 45mm/43mm. Packing: pouch – 50pcs., carton 2000 pcs.
FES-X4543B		Adult, monitoring, carbon snap, Ag /AgCl sensor, RTG, TM, MR conditional, solid gel, PE foam, size 45mm/43mm. Packing: pouch – 50pcs., carton 2000 pcs.
FES-55A		Adult, Holter, Stress/monitoring, metal snap, Ag /AgCl sensor , solid gel, PE foam, dia 55mm. Packing: pouch – 50pcs., carton 2000 pcs.
FES-50A		Adult, Holter, Stress/monitoring, metal snap, Ag/AgCl sensor, solid gel, PE foam, dia 50mm. Packing: pouch – 50pcs., carton 2000 pcs.

FES-X50A		Adult, Holter, Stress/monitoring, carbon snap, Ag/AgCl sensor, RTG, TM, MR conditional, solid gel, PE foam, dia 50mm. Packing: pouch – 50pcs., carton 2000 pcs.
FEL-55A		Adult, Holter, Stress/monitoring, metal snap, Ag/AgCl sensor, wet gel, PE foam, dia 55mm. Packing: pouch – 50pcs., carton 1500 pcs.
FEL-50A		Adult, Holter, Stress/monitoring, metal snap, Ag/AgCl sensor, wet gel, PE foam, dia 50mm. Packing: pouch – 50pcs., carton 1500 pcs.
CES-55A		Adult, Holter / Stress, monitoring, metal snap, Ag/AgCl sensor, solid gel, nonwoven, dia. 55mm, Packing: pouch – 50pcs., carton 2000 pcs.
CES-50A		Adult, Holter / Stress, monitoring, metal snap, Ag/AgCl sensor, solid gel, nonwoven, dia. 50mm. Packing: pouch – 50pcs., carton 2000 pcs.
FES-XB5136C		Adult, monitoring, ø 4mm connector (BANANA), Ag/AgCl sensor, X-Ray, CT, MR conditional, solid gel, PE foam, size 51/36 mm. Packing: pouch – 30 pcs., carton 900 pcs.
FES-4137B		Pediatric, Holter/monitoring, metal snap, Ag/AgCl sensor, PE foam, solid gel, size 41mm x 37mm. Packing: pouch – 50pcs., carton 2500 pcs.
FES-3631B		Pediatric/New born, monitoring, metal snap, Ag/AgCl sensor, PE foam, solid gel, size 36 mm x 31 mm. Packing: pouch – 50pcs., carton 2500 pcs.

<p>FES-30A</p>		<p>Pediatric/New born, monitoring, metal snap, Ag/AgCl sensor, PE foam, solid gel, dia. 30mm. Packing: pouch – 50pcs., carton 2500 pcs.</p>
<p>FES-25A</p>		<p>New born, monitoring, metal snap, Ag/AgCl sensor, PE foam, solid gel, dia. 25mm. Packing: pouch – 50pcs., carton 2500 pcs.</p>

FARUM sp. z o.o.
74, Jagiellonska Street, 03-301 Warsaw, Poland

Contact data:

Anna Pustoła, Export & Import Executive

Phone/Fax: (+48 22) 811-19-22

E-mail: export@farum.com.pl

Web: www.farum.pl