





## AndraBalloon<sup>®</sup> PTA



- Optimized for Stent dilatation
- Balloon diameter from 10 up to 32 mm
- Balloon length from 20 up to 60 mm

# AndraBalloon® PTA-Catheter



- Semi-compliant balloon made of PA
- Optimized folding for better rewrapping and easy repassage
- Coaxial shaft design for fast inflation and deflation
- Optimized PTA-Balloon for Stent dilatation
- 2 Pt/Ir X-Ray markers on the inner shaft

## Order Codes AndraBalloon® PTA-Catheter

Order Code	Balloon Diameter (mm)	Balloon Length (mm)	Nominal Pressure	Rated Burst Pressure	Introducer Sheath
ABL 10/20-1200 ABL 10/30-1200 ABL 10/40-1200	10	20/30/40	400kPa (4bar)	1000kPa (10bar)	7F
ABL 12/20-1200 ABL 12/30-1200 ABL 12/40-1200	12	20/30/40	400kPa (4bar)	1000kPa (10bar)	7F
ABXL 14/30-1200 ABXL 14/40-1200 ABXL 14/60-1200	14	30/40/60	400kPa (4bar)	800kPa (8bar)	8F
ABXL 16/30-1200 ABXL 16/40-1200 ABXL 16/60-1200	16	30/40/60	400kPa (4bar)	800kPa (8bar)	10F
ABXL 18/20-1200 ABXL 18/30-1200 ABXL 18/40-1200 ABXL 18/60-1200	18	20/30/40/60	400kPa (4bar)	700kPa (7bar)	10F
ABXL 20/30-1200 ABXL 20/40-1200 ABXL 20/60-1200	20	30/40/60	400kPa (4bar)	700kPa (7bar)	11F
ABXL 22/40-1200 ABXL 22/60-1200	22	40/60	400kPa (4bar)	700kPa (7bar)	11F
ABXL 24/40-1200 ABXL 24/60-1200	24	40/60	400kPa (4bar)	700kPa (7bar)	12F
ABXL 26/40-1200 ABXL 26/60-1200	26	40/60	400kPa (4bar)	500kPa (5bar)	13F
ABXL 28/40-1200 ABXL 28/60-1200	28	40/60	400kPa (4bar)	500kPa (5bar)	13F
ABXL 30/40-1200 ABXL 30/60-1200	30	40/60	400kPa (4bar)	500kPa (5bar)	13F
ABXL 32/40-1200 ABXL 32/60-1200	32	40/60	400kPa (4bar)	500kPa (5bar)	13F

