





## AndraSnare<sup>®</sup> & AndraSnare Micro<sup>®</sup>



- Gold plated atraumatic 3D Nitinol Snare
- Super-elastic & self-expanding
- 1:1 control of the loop

# AndraSnare® & AndraSnare Micro®



- 1 Atraumatic Nitinol loop with gold coating for better visualization in the X-rays
- 2 Pre-formed tip for better navigation
- 3 X-ray marker at the shaft output
- 4 Shaft visible in X-rays
- 5 Bend protection for shaft reinforcement at the Luer
- 6 Luer connection
- 7 Insertion part for simple insertion of the loop into the catheter
- 8 Large torquer for good operation
- 9 Continuous Nitinol wire for 1:1 control of the loop

## Order Information AndraSnare®

Product Code	Loop- Diameter	Length (cm) Loop/Guidewire	Introducer Sheath (cm)	Introducer Sheath (F)
AS-5	5 mm (3-5 mm vessel)	125	110	4
AS-10	10 mm (6-10 mm vessel)	125	110	4
AS-15	15 mm (11-15 mm vessel)	125	110	4
AS-20	20 mm (16-20 mm vessel)	125	110	4
AS-25	25 mm (21-25 mm vessel)	125	110	5
AS-25/65	25 mm (21-25 mm vessel)	65	50	5
AS-25/150	25 mm (21-25 mm vessel)	150	135	5
AS-30	30 mm (26-30 mm vessel)	125	110	5
AS-35	35 mm (31-35 mm vessel)	125	110	5

## Order Information AndraSnare® for small vessels

Product Code	Loop- Diameter	Length (cm) Loop/Guidewire	Introducer Sheath (cm)	Introducer Sheath (F)
ASM-2	2 mm (1-2 mm vessel)	175	150	2,3 - 3
ASM-4	4 mm (3-4 mm vessel)	175	150	2,3 - 3
ASM-7	7 mm (5-7 mm vessel)	175	150	2,3 - 3

