





AndraBasket[®]




- Single-handed operation
- Removable handgrip / sheath
- Torqueable 1:1
- Excellent pushability

AndraBasket®



- 1 Atraumatic tip for safty placement of the basket
- 2 Helical nitinol basket create ideal catching conditions
- 3 Flexible shaft
- 4 Bend protection at the Luer
- 5 Handle permits simple application – One hand technique
- 6 Luer Port Connector



Order Information AndraBasket® and AndraBasket Micro®

Product Code	French Size	Basket Diameter	Number of Struts	Sheath Lenght (cm)
ABT-H-4-2-12-090H	2	12 mm	4	90
ABT-H-4-2.5-12-090H	2,5	12 mm	4	90
ABT-H-4-2.5-12-120H	2,5	12 mm	4	120
ABT-H-4-3-15-090H	3	15 mm	4	90
ABT-H-4-3-15-120H	3	15 mm	4	120
ABT-H-4-4-17-090H	4	17 mm	4	90
ABT-H-4-4-17-120H	4	17 mm	4	120

Order Information AndraBasket Gastro®

Product Code	French Size	Basket Diameter	Number of Struts	Sheath Lenght (cm)
ABT-G-4-5-25-215H	5	25 mm	4	215
ABT-G-4-5-30-215H	5	30 mm	4	215
ABT-G-4-7-25-215H	7	25 mm	4	215
ABT-G-4-7-30-215H	7	30 mm	4	215

