

AndraStent XL & XXL Balloon Expandable Stent System




- Cobalt Chromium Stent - Performing High Performance & Minimal Puncture Site
- "Hybrid Cell Design" For High Flexibility
- Very High Product Variability



DELIVERING LESS INVASIVE FUTURE TECHNOLOGY

AndraStent XL & AndraStent XXL



- The AndraStent XL & the AndraStent XXL have been specifically designed to meet the high requirements for extra large vessels & procedures where XL & XXL Stents are needed.
- Improved flexibility to the conventional XL Stents allow much better delivery performance
- New materials deliver lower crimping profiles with lower sheath compatibilities but remains high radial strength

Technical Specifications:

Material:	HD Cobalt Chromium MP35N
Stent Design:	Hybrid Cell Design
Dilatation Range Andramed XL-Stent:	15mm – 25mm Ø
Recommended Sheath XL-Stent:	8 F (+ 1F) depending on the PTA
Dilatation Range Andramed XXL Stent:	20mm – 32mm Ø
Recommended Sheath XXL Stent:	10F (+ 1F) depending on the PTA
Stent Recoil:	< 2.8%
Surface / Artery Ratio	18% -23%
Strut-thickness:	0.010"

Foreshortening Table:

Stent length (mm)	Expansion Diameter (ø mm)					Product Code
	Expanded Length AndraStent XL (mm)					
	ø 15	ø 18	ø 20	ø 22	ø 25	
13	11	10	10	8	7	AS13XL
17	15	13	13	11	10	AS17XL
21	19	17	16	14	13	AS21XL
26	22	20	20	17	15	AS26XL
30	26	24	23	20	18	AS30XL
35	30	27	26	23	21	AS35XL
39	34	31	30	26	24	AS39XL
43	38	34	33	29	26	AS43XL
48	41	38	36	32	29	AS48XL
57	49	45	43	38	34	AS57XL

Stent length (mm)	Expansion Diameter (ø mm)					Product Code
	Expanded Length AndraStent XXL (mm)					
	ø 20	ø 22	ø 25	ø 28	ø 32	
17	16	15	14	12	10	AS17XXL
21	20	19	18	15	13	AS21XXL
26	24	22	21	18	16	AS26XXL
30	28	26	25	22	18	AS30XXL
35	33	30	28	25	21	AS35XXL
39	36	34	32	28	24	AS39XXL
43	41	38	36	31	26	AS43XXL
48	45	42	39	34	29	AS48XXL
57	53	49	46	40	34	AS57XXL

Order Informations:

Order Codes AndraStent XL and XXL										
Stent Length										
	13 mm	17 mm	21 mm	26 mm	30 mm	35mm	39 mm	43 mm	48 mm	57 mm
Andrastent XL	AS13XL	AS17XL	AS21XL	AS26XL	AS30XL	AS35XL	AS39XL	AS43XL	AS48XL	AS57XL
Andrastent XXL	-----	AS17XXL	AS21XXL	AS26XXL	AS30XXL	AS35XXL	AS39XXL	AS43XXL	AS48XXL	AS57XXL

Included: Metal Introducer A-INT1213F



Your Local Distributor:



Andramed GmbH
 Silvanerstr. 8
 D-72766 Reutlingen
 ph: +49(0)7127 972290
 Fax: +49(0)7127 972291
 www.andramed.com

Rev. 6 08.08.2010